CEL✚X[™] PPH

INSTRUCTIONS FOR USE

en

Uterine Haemostatic Tamponade

- INTENDED PURPOSE:**
CELOX[™] PPH is intended to be a physical haemostat for control of emergency bleeding.
- INDICATION FOR USE:**
CELOX[™] PPH is intended to be applied and removed vaginally to provide control and treatment of uterine postpartum haemorrhage when conservative management is warranted.
- DESCRIPTION:**
CELOX[™] PPH is a uterine haemostatic tamponade provided in a z-fold/concertina format. It has CELOX[™] haemostatic granules loosely adhered to a fabric (3m x 7.6cm).
- PATIENT GROUP:**
Females of childbearing age.
- USERS:**
The product is intended for use by physicians and midwives trained and experienced in obstetrical techniques.
- CLINICAL BENEFIT:**
Postpartum haemorrhage can be life-threatening. The use of the device will reduce the likelihood of needing hysterectomy (or other interventional methods) by controlling postpartum haemorrhage.
- CONTRAINDICATIONS**
- Subjects who present with uterine rupture.
 - Unresolved uterine inversion.
 - Current cervical cancer.
 - Current purulent infection of the vagina, cervix, uterus.
 - Ongoing Pregnancy.
- WARNINGS**
- Do not use if the primary pack is damaged or open as the device may not be sterile.
 - Keep away from liquids prior to use. The CELOX[™] granules will gel in contact with liquid, making application difficult and reduce haemostatic effectiveness.
 - As with other haemostatic agents, do not aspirate blood into cell saver equipment or autologous blood salvage circuits. There is a potential risk of embolism if the blood is returned to the patient.
 - The haemostatic effect of CELOX[™] PPH is not enhanced by the addition of topical thrombin, the activity of which may be destroyed by the pH of CELOX[™] PPH.
 - CELOX[™] PPH must not be used as a vaginal/cervical pack for uterine postpartum haemorrhage as uncontrolled bleeding could continue in the uterus.
 - For single use only. Re-use could potentially result in the risk of infection and lack of device performance.
 - If the patient develops high fever (>38°C) following haemostasis remove the device.
 - Do not leave CELOX[™] PPH for longer than 24h in the body due to possible risk of infection.
 - Transvaginal application of CELOX[™] PPH may not be possible when cervix dilation has not started e.g., in case of elective c-section.
 - Do not use suturing techniques (e.g., uterotomy closure, B-Lynch, Pereira) whilst CELOX[™] PPH is in the uterus as this may result in CELOX[™] PPH being incorporated into the suture causing difficulty for removal and risk of product tearing.
- PRECAUTIONS**
- If used concurrently with compression suturing e.g., B-Lynch, the sutures should be applied prior to CELOX[™] PPH application to avoid excessive compression by the procedure resulting in difficulty to remove the CELOX[™] PPH after treatment completion.
 - Caution is advised when suturing (e.g., for vaginal or cervical tears) that CELOX[™] PPH is not included in the stitch, which upon removal of CELOX[™] PPH may tear leaving a fragment attached to the stitch.
 - During use avoid excessive manipulation of the patient when transferring to hospital or between departments as this may disturb the device, potentially resulting in rebled.
- STORAGE**
There are no special storage conditions required for the device, except to keep the shipper cases and instruction leaflet dry.
- INSTRUCTIONS FOR USE**
- CELOX[™] PPH should be used when postpartum haemorrhage does not respond to initial management such as uterine massage, volume replacement, and standard medical treatment with uterotonic and antifibrinolytic agents.
 - Tear open the package at the tear notches. Sterile technique should be observed in delivering the sterile CELOX[™] PPH to the field of application.
 - After caesarian deliveries ensure the uterotomy is closed prior to applying the device transvaginally.
 - Ensure any retained placenta tissue, blood clots or membranes are removed from the uterus prior to uterine packing with CELOX[™] PPH.
 - Remove any equipment applying external pressure to the uterus (e.g., pressure pads) prior to applying the device.
 - Use at least one CELOX PPH for each patient and do not cut or tear CELOX PPH prior to application as this could increase the risk of material fragments being left in the uterus. Use a maximum of 2 CELOX PPH per patient over the course of treatment.
 - Allow the CELOX PPH to unfold externally to the patient whilst applying the device.
 - Place the patient in a lithotomy position, and with the aid of a speculum grasp the cervix with an atraumatic instrument, if necessary, and manually insert the unfolding CELOX[™] PPH into the uterine cavity with traumatic forceps up to the fundus. Avoid excessive force and use ultrasound guidance if available to guide application to ensure the device reaches the site(s) of bleeding and to minimise the risk of damaging the uterine wall. If ultrasound is not available, use external manual palpation to avoid injury to the uterine wall by the insertion device.
 - CELOX[™] PPH must contact the site(s) of bleeding.
 - CELOX[™] PPH should be used to cover the uterine wall and if needed use secondary sterile dressings/gauze (non CELOX[™]) for packing. Leave sufficient length of the end of each CELOX[™] PPH protruding through the cervix and vagina to allow subsequent removal.
 - If needed, the vagina may be packed with secondary sterile dressings/gauze (non CELOX[™]) for compression of ascending vessels.
 - While CELOX[™] PPH is in place, compression is to be applied for a sufficient duration according to clinical judgement to enable haemostasis.
 - If postpartum haemorrhage continues, remove CELOX[™] PPH and revert to the next stage of the standard of care/protocol for post-partum haemorrhage.
 - Following haemostasis, excess CELOX[™] PPH can be cut from the ends protruding from the vagina (leaving sufficient length of each CELOX[™] PPH to be noticeable and to allow removal).
 - Dispose of any excess material and packaging as per local clinical waste management guidelines.
 - Use local standards of care to identify the number of CELOX[™] PPH that have been applied

- to the patient e.g., use of wristbands. Identify in the patient file that CELOX[™] PPH has been used and indicate when it needs to be removed. Ensure the patient is informed CELOX[™] PPH is in place and needs to be removed by a professional within 24hrs.
- Monitor the patient carefully for signs of rebleeding (e.g. clinical signs, monitoring of blood values and using ultrasound scans if available). Treatment with oxytocic/uterotonic drugs may continue during this time.
 - Closely monitor for signs disseminated intravascular coagulation (DIC), and in such cases, emergency intervention per hospital protocol should be followed.
 - Leave CELOX[™] PPH in place for 24 hours if clinically feasible, but no longer.
 - CELOX[™] PPH should be removable by pulling the end left outside the vagina.
 - Upon removal, inspect the end of the fabric for an intact edge to determine integrity and completeness. Adequate revision action steps should be initiated when incomplete removal is observed or suspected. No fabric fragments should be left behind.
 - Dispose of the used CELOX[™] PPH according to standard hospital procedures for clinical waste.
 - It is expected loose residual granules or gelled material from CELOX[™] PPH will be expelled in the lochia of the patient and potentially upon menstruation restarting. This could potentially be several weeks after use. This information should be conveyed to the patient.

- ADDITIONAL INFORMATION**
- Not made with natural rubber latex.
 - The device is effective in patients that are coagulopathic.
 - CELOX[™] PPH and residuals are not intended to be absorbed by the body.
 - The safety and effectiveness of CELOX[™] PPH with other topical haemostatic agents or sealants/adhesives has not been established.
 - Contains chitosan from shellfish, allergy studies showed no adverse reactions to similar chitosan.
 - The Summary of Safety and Clinical Performance (SSCP) for the device is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. The Eudamed website is <https://ec.europa.eu/tools/eudamed> and the Basic UDI-DI is 506020663088381712L.
 - Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Medtrade Products Ltd, Electra House,
Crewe Business Park, Crewe, CW1 6GL, UK.
www.celoxpph.com

MedEnvoy, Prinses Margrietplantsoen
33 – Suite 123, 2595 AM The Hague, The Netherlands

FG08838321
Artwork Ref: MTP-24-2767 Rev 1
Issue Date: 15/Jan/2025
CELOX[™] is a trademark of Medtrade Products Ltd. All rights reserved.

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